



System Zeramex XT

 <p>RB16550</p>  <p>XT Abutment RB Straight</p>  <p>XT Abutment WB Straight</p>  <p>XT Abutment RB Angular</p>  <p>XT Abutment WB Angular</p> <p>Anzugsdrehmoment: 25Ncm</p>	 <p>RB36553</p>  <p>XT Scanbody RB</p>  <p>XT Scanbody WB</p>  <p>Alle XT Abutments RB & WB</p> <p>Anzugsdrehmoment: 5Ncm</p>	 <p>RB36550</p>  <p>XT Provisorium RB</p>  <p>XT Provisorium WB</p> <p>Anzugsdrehmoment: 15Ncm</p>
---	--	---

System Zeramex P6

 <p>P15450</p>  <p>P6 Abutment SN Straight</p>  <p>P6 ZERABASE SN</p> <p>Anzugsdrehmoment: 25Ncm</p>	 <p>P15451</p>  <p>P6 Abutment SN Angular</p> <p>Anzugsdrehmoment: 25Ncm</p>	 <p>P35550</p>  <p>P6 Provisorium SN</p>  <p>P6 Scanbody SN</p> <p>Anzugsdrehmoment: 15Ncm</p>
 <p>P16450</p>  <p>P6 Abutment RN Straight</p>  <p>P6 ZERABASE RN</p> <p>Anzugsdrehmoment: 35Ncm</p>	 <p>P16451</p>  <p>P6 Abutment RN Angular</p> <p>Anzugsdrehmoment: 35Ncm</p>	 <p>P36550</p>  <p>P6 Provisorium RN</p>  <p>P6 Scanbody RN</p> <p>Anzugsdrehmoment: 15Ncm</p>